

Product Name:	PCR Tube			
Product Application:	NEST has both transparent and white PCR strip tubes, which respectively apply to the general PCR and qPCR reactions.			
Product Materials:	polypropylene (PP).			
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 5 years after date of production (ensure package is in good)			
Sterilization:	Yes (some of all catalog numbers). Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137			
Pyrogens:	Non-Pyrogenic			
RNase/DNase:	DNase/RNase free			
BSE/TSE:	These products are deemed animal free.			
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Vacuum Seal test: qualified Freeze testing: qualified Centrifugation test: qualified Evaporation rate test: qualified Autoclave test : qualified Packaging inspection: qualified			
Features:	1. Made from USP grade VI polypropylene 2. 0.1mL tubes are low profile which are designed to minimize the condensate contamination, particularly for the overnight experiment, and also reduce the evaporation of the reaction solution. 3. Universally applicable to the mainstream PCR instrument and real-time qPCR on the market. Thin-walled design produces high thermal conductivity, enabling the reaction solution in the tube to reach the target temperature as quickly as possible 4. The tube cap has excellent sealing performance, and the capping is smooth, so that the loss of reaction volume is less than 5% when using a heated lid 5. Free of Dnase and Rnase, human DNA and PCR inhibitor. 6. Autoclavable (121℃, 20min) 7. White strip tubes increase the fluorescent signal strength of real-time PCR.			

Product Range

Cat. No.	Product Description	Sterilization	/Pack	/Case
401001	0.2 mL PCR Tubes, Flat Caps, Clear	No	500	10000
403102	0.1 mL 8-well PCR Tube Strips, Clear	No	125	1250
403122	0.1 mL 8-well PCR Tube Strips, Clear	Yes	125	1250
403002	0.2 mL 8-well PCR Tube Strips, Clear	No	125	1250
403022	0.2 mL 8-well PCR Tube Strips, Clear	Yes	125	1250

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea


NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137


403112	0.1 mL 8-well PCR Tube Strips, White	No	125	1250
403132	0.1 mL 8-well PCR Tube Strips, White	Yes	125	1250
403012	0.2 mL 8-well PCR Tube Strips, White	No	125	1250
403032	0.2 mL 8-well PCR Tube Strips, White	Yes	125	1250
404001	0.2 mL 8-well PCR Tube Strips, With Flat caps, Clear	No	12	1200

Technical Drawing

401001



403102/403112




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

403002/403012



404001



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)