

Product Name:	Petri Dish
Product Application:	Suitable for bacterial culture.
Product Materials:	Polystyrene (PS)
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Packaging inspection: qualified
Features	1. High clarity polystyrene 2. Flat transparent surface 3. Numeric mark on the bottom for easy identification 4. PA/PE packaging 5. Vent for gas exchange 6. Stackable for easy storage and handling 7. Non-pyrogenic 8. Sterilized by E-Beam, SAL=10⁻⁶

Product Range

Cat. No.	Spec(mm)		Volume (mL)	/Pack	/Case
	Diameter	Height			
706011	35	12	5	20	500
754001	60	15	15	20	500
722011	65	15	13-15 (RODAC Dish)	20	500
752001	90	15	40	20	500
752002 (Double Packaged)	90	15	40	20	500
752003	90	15	40	5	500
752004	90	15	40	10	500
752011	90	15	20 x 2 compartments	20	500
752021	90	15	13 x 3 compartments	20	500
752031	90	15	10 x 4 compartments	20	500
753401	90	20	40	20	300

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

752501	90	25	40	13	325
715011	150	15	60	10	100
755001	150	25	130	5	100
753901	100 x100 x15 (Square)	17.25	45	20	500


Technical Drawing

35 mm Petri Dishes




Cap


60 mm Petri Dishes



Cap



Bottom



Bottom


Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

65 mm RODAC Dish



90*15 mm Petri Dishe




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)


90mm Petri Dishes(2 compartments)



90mm Petri Dishes(3 compartments)



90mm Petri Dishes(4 compartments)




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

90*20 mm Petri Dishes



90*25 mm Petri Dishes




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

100*15 mm Petri Dishes



150*15 mm Petri Dishes




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

150*25 mm Petri Dishes



100*100*15 mm Square Petri Dishes



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)