


Product Name:	L-Spreader
Product Application:	Applied in microbiology experiments. For inoculating bacteria on agar medium.
Product Materials:	Acrylonitrile Butadiene Styrene (ABS)
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens: RNase/DNase	Non-Pyrogenic DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Packaging inspection: qualified
Features	1. Smooth and flat surface that avoids problems such as contamination or uneven application of spread plate in bio-labs. 2. Length of handle: 144mm, diagram of spreading: 33mm, suitable for 3.5cm, 6cm, 9cm, 10cm petri dishes, etc. Ideal for spreading cells or bacteria evenly in the culture dish. 3. Easy-tear sterile individual packaging.

Product Range

Cat. No.	Product Description	/Pack	/Case	Pcs/Case
711001	L-Spreader, Sterile	200	10	2000
711002	L-Spreader, Sterile	5	100	500

Technical Drawing



Head Office Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)