

Product Name:	Cryo Boxes
Product Application:	For centralized placement of cryogenic vials.
Product Materials:	PC (cap) + ABS (bottom)
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)
Sterilization:	No
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance Inspection: qualified Packaging Inspection: qualified
Features	1. The cap is made of polycarbonate plastic (PC), and the bottom is made of acrylonitrile butadiene styrene plastic (ABS). 2. Clear polycarbonate cover for easy identification. 3. Individually packed. 4. Alphanumeric grid for easy sample identification. 5. Each box includes drain holes and air vents to minimize condensation. 6. Withstand temperature range: Approx.-196°C to 70°C.

Product Range

Cryogenic Box			
Cat. No.	Product Description	/Pack	/Case
For 2D Cryogenic Tubes			
616651	For 1.5 mL vials	10*10	14
616051	For 2.0 mL vials/1.8 mL Internal thread vials	10*10	12
616151	For 4.0 mL vials	10*10	8
616551	For 5.0 mL vials	10*10	6
For SBS Format Cryogenic Tubes			
616211	For 0.5 mL external thread vials	12*8	10
616212	For 0.5 mL internal thread vials	12*8	10
616711	For 0.75 mL vials	12*8	10
616361	For 1.0 mL external thread vials	12*8	10
616362	For 1.0 mL Internal thread vials	12*8	10

Head Office Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

616212	For 1.5 mL vials	12*8	10
616041	For 2.0 mL vials	8*6	10
616141	For 4.0 mL vials	8*6	10
616841	For 6.0 mL vials	6*4	10
616891	For 8.0 mL vials	6*4	10

Cap Inserts

Cat. No.	Color	Sterilization	/Pack	/Case
611201	White	No	100	1000
611202	Red	No	100	1000
611203	Yellow	No	100	1000
611204	Blue	No	100	1000
611205	Green	No	100	1000
611206	Purple	No	100	1000


Note:

New Cap Inserts are compatible with new cryogenic tubes.

New Cap Inserts do not match 1mL, 0.75mL, 0.5mL cryogenic tubes, but all other sizes do.

Technical Drawing

Cap Inserts






Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Technical Drawing



10*10 Cryo boxes for 1.5 mL vials







10*10 Cryo boxes for 2.0/1.8 mL vials

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com





Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)



10*10 Cryo boxes for 4.0 mL vials

10*10 Cryo boxes for 5.0 mL vials



12*8 Cryo boxes for 0.5 mL external thread vials



12*8 Cryo boxes for 0.5 mL internal thread vials

Head Office



Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea



NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)





12*8 Cryo boxes for 0.75 mL vials



8*6 Cryo boxes for 2.0 mL vials



12*8 Cryo boxes for 1.0 mL external thread vials





12*8 Cryo boxes for 1.0 mL internal thread vials

Head Office



Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea



NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)



8*6 Cryo boxes for 4.0 mL vials



6*4 Cryo boxes for 6.0 mL vials



6*4 Cryo boxes for 8.0 mL vials

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)