

Supelco.

1.14660.0001

50 Tests

MQuant®

Iron Test

0.05 - 1.0 mg/l Fe

Merck KGaA, 64271 Darmstadt, Germany, Tel. +49(0)6151 72-2440
www.sigmaaldrich.com

EMD Millipore Corporation, Burlington MA 01803, USA,
Tel. +1-978-715-4321

MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
Ontario, L6H 6J8, Canada, Phone: +1 800-565-1400

7.69146.6006-XXXXXX 06/2024

5 ml

1 x Fe-1

4 x Fe-2

3 min

0.4

0.6

0.8

1.0

mg/l Fe

0

0.05

0.1

0.2

mg/l Fe