



Reference: 064-BA7069

Scharlau Microbiology - Technical Data

Product: PBS Buff pH= 7,2

Specification

Phosphate solution used in microbiological procedures.

Presentation

10 Prepared bottle
Bottle 500 ml
with: 300 ± 5 ml

Packaging Details

1 box with 10 bottles 500 ml. Injectable cap: Plastic screw inner cap. The use of syringes needles with a diameter greater than 0.8 mm is not recommended.

Shelf Life

12 months

Storage

8-25 °C

Composition

Composition (g/l):

Potassium dihydrogen phosphate..... 0.20
Disodium hydrogen phosphate..... 1.15
Sodium chloride..... 8.00
Potassium chloride..... 0.20

Description /Technique

Buffered diluent used in routine processes of preparation of dilution banks for several microbiological tests.

Inoculate according to final purpose, samples and validated methods.

Quality control

Physical/Chemical control

Color : Colourless

Microbiological control

According to USP & European Pharmacopeia

Inoculate with 10-100 CFU according to harmonized Pharmacopoeiae or with 10^4 - 10^6 CFU for selectivity.

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Subculture onto appropriate culture media after holding at 20-25 °C for 45 min. to 1 h.

Aerobic. Incubation at 30-35 °C for 18-72h (bacteria) and 20-25 °C for 3-5 days (moulds and yeast).

Microorganism

Aspergillus brasiliensis ATCC® 16404, WDCM 00053

Bacillus subtilis ATCC® 6633, WDCM 00003

Candida albicans ATCC® 10231, WDCM 00054

Escherichia coli ATCC® 8739, WDCM 00012

Ps. aeruginosa ATCC® 9027, WDCM 00026

Staphylococcus aureus ATCC® 6538, WDCM 00032

Growth

Good. Recovery ±30% T0 (original enumeration)

Good. Recovery ±30% T0 (original enumeration)

Good. Recovery ±30% T0 (original enumeration)

Good. Recovery ±30% T0 (original enumeration)

Good. Recovery ±30% T0 (original enumeration)

Good. Recovery ±30% T0 (original enumeration)

Sterility control

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Bibliography

- ATLAS, R.M. & L.C. PARKS (1997) Handbook of microbiological media. CRC Press. BocaRaton .Fla. USA.
- Dulbecco and Vogt (1954) J. Exp. Med. 99. 167-182.
- Paul J. (1965) 'Cell and Tissue Culture' 3rd ed, Livingstone Ltd., London.