

Product Specifications

CL02.1443

Sodium hydroxide 32 weight % solution a.r.

Sódio hidróxido 32 p/p % solução p.a.
Sodio idrossido 32 p/p % soluzione p.a.
Sodio hidróxido 32 p/p % solución p.a.
Natronlauge 32 w/w % Lösung z..
Sodium hydroxyde 32 poids % solution p.a.
Natriumhydroxide 32 gew. % oplossing p.a.
Kjeldahl N Det., Ph. Eur.
432 g NaOH / l H₂O

Mol weight: 40.00 g/mol
Density: 1,349 g/ml
HS Nr: 28151200
CAS Nr: 1310-73-2
UN: 1824
ADR: 8,II



Specifications:

Assay	: 31 - 33%
Sodium Carbonate	<0.5%
Aluminium (Al)	<0.0005%
Iron (Fe)	<0.0005%
Heavy Metals as Lead (Pb)	<0.0002%
Chloride (Cl)	<0.001%
Nitrogen (N)	<0.0001%
Sulfate (SO ₄)	<0.005%
Sulfur (S)	<0.0001%
Cyanide (CN)	<0.0001%
Silicate (SiO ₂)	<0.01%