

Technical Data Sheet

Product information

Product code	B601.98.0005
Product name	Buffered Peptone Water
Short name	BPW
Packing description (size, type, filling, amount)	15mm, Swabs, 5ml, 25pc(s)

Intended user	For professional use only.
---------------	----------------------------

Storage Conditions

Storage temperature	2-25°C
Storage conditions	Keep the products in a low humidity environment, protected from moisture and light.
Shelf Life	See expiry date on the label (valid for closed packaging).
Shelf life after production	26 Weeks

Typical Composition *)

Ingredient	Amount of	Per liter
Enzymatic digest of casein	10	g
Sodium chloride	5	g
Di-sodium hydrogen phosphate	3,5	g
Potassium dihydrogenphosphate	1,5	g
Demiwater	1000	ml

*) Adjusted as required to meet performance standards.

Analytical Performance characteristics

pH	7,0 +/- 0,2
Colour	yellow
Clear	Yes
Particles	No

Excellent
Environment
& Life

ONE¹
XDG

Microbiological Performance

Conditions:

Name	Incubation time	Incubation temp (°C)	Environment	Transfer	Incubation time	Incubation temp (°C)	Environment
inc-A	16 to 20 hours	37 +/-1	Aerobic				

Name	Microorganism	Identification number	Method	Specification	Electivity
inc-A	Salmonella Enteritidis	ATCC 13076 WDCM 00030	Productivity of liquid media	Slight to good growth	
inc-A	Salmonella Typhimurium	ATCC 14028 WDCM 00031 NCTC 12023	Productivity of liquid media	Slight to good growth	

Microbiological Contamination

	Incubation time	Incubation temp (°C)	Transfer	Incubation time	Incubation temp (°C)	Specification
Inc-a	44 to 52 hours	30 +/- 1				No growth

Excellent
Environment
& Life

ONE¹
XDG

Contact:

Xebios Diagnostics Group GmbH
Tellerlingstrasse 49
45097 Düsseldorf, Germany
www.xebios.com
info@xebios.com
tel: +49 211 78 88 07 - 00

Date of creating: 08-08-2024

Version number: 2



Excellent
Environment
& Life

ONE¹
XDG