

Solstice

Honeywell
Specialty Chemicals
Seelze GmbH
P.O. Box 10 02 62
30918 Seelze /Germany
Tel.: 0049 (0)5137-999-0
Fax.: 0049 (0)5137-999-170

SPECIFICATION: 03652 1-Octanol CHROMASOLV, for HPLC, 99%

Print Date: 02.02.2026 / Revision index: 01 / Effective: 07.10.2016

	complying		
identity (IR)			
assay (GC)	min.	99.0	%
non-volatile matter	max.	0.004	%
water (Karl Fischer)	max.	0.01	%
absorbance at 400 nm	max.	0.01	
absorbance at 300 nm	max.	0.01	
absorbance at 250 nm	max.	0.05	
absorbance at 235 nm	max.	0.20	
absorbance at 225 nm	max.	0.50	
absorbance at 215 nm	max.	1.0	
APHA	max.	10	

The above specification shows our present quality of the material.

We reserve the right to change this specification, if necessary.

We warrant the quality of our products according to our General Conditions of Sale.