



Acide tannique

Synonyms:

Gallotannin

C76H52O46

Classification

Danger

Molecular weight

1701.2
3 ONU

CAS

1401-55-4 Transport hazard
class:

EEC-N:

215-753-2 Packing

Acide tannique

>

ERBApharm - Selon pharmacopée: Ph.Eur.-
USP-FU

Grade

Description	Yellowish powder	-
Identification (I.R.)	Positive	-
Appearance of solution	Conform*	Ph.Eur.
Loss on drying	<= 12.0	%
Sulfated ashes	<= 0.1*	%
Dextrin, gum, salt, sugar	Conform*	Ph.Eur.



Acide tannique > ERBAPharm - Selon pharmacopée: Ph.Eur.-USP-FU		Grade
Resin	Conform*	Ph.Eur.
Heavy metals (as Pb)	<= 40*	ppm

ID	Size	Packaging
0000000000000307152	5 kg	Seau plastique
0000000000000307153	25 kg	Seau plastique
0000000000000307157	1 kg	Sachet



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val de
 Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
 operates with a Certified
 Quality Management
 System

