



Ethanol absolu anhydre

Synonyms:

C₂H₅OH**Classification****Danger**

Molecular weight

46 ONU**1170**

CAS

64-17-5 Transport hazard
class:**3**

EEC-N:

200-578-6 Packing**II****Ethanol absolu anhydre** > **RS - Pour histologie - CE - IVD****Grade Specificque**

Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Density at 20°C	0.790 - 0.793	-
Water (K.F.)	<= 0.1	%
Residue on evaporation	<= 30	ppm
Assay (alcoholic) at 20°C	>= 99.9	% (V/V)
Description	Clear colourless liquid	-



Identification (I.R.)	Positive	-
Density at 20°C	0.790 - 0.793	-
Water (K.F.)	<= 0.1	%
Residue on evaporation	<= 30	ppm
Assay (alcoholic) at 20°C	>= 99.9	% (V/V)

ID	Size	Packaging
000000000000308715	5 l	Bidon plastique