



Product Data Sheet

LLG-PE Snap Ring Seals ND11, ready assembled

Cat. No.: 4.001 544



Article description: LLG-11 mm PE Snap Ring Cap, transparent, centre hole, soft Version. Silicone white/PTFE red, 45°, shore A, 1,0 mm

Cap size: ND11

Colour: Transparent

Material septa: Silicone/PTFE

Colour septa: White/red

Shore hardness: 45 A

Slit: Without

Packaging: pack of 100 pcs

This document was sent online and is valid without signature