



Product data sheet

LLG-Aluminium crimp caps N 8

Art.-Nr.: 4.008 198

Cap

Typ of Thread	Crimp Cap Aluminium
Material	clear lacquered
Colour	

Dimensions

Skirt Height in mm	$4,95 \pm 0,1$
Material Thickness in mm	$0,20 \pm 0,01$
Inner Diameter in mm	$7,80 \pm 0,05$
Outer Diameter in mm	$8,20 \pm 0,07$
Size of Centre Hole in mm	$4,00 \pm 0,2$

Septum

Material and colour	PTFE red / Silicon white / PTFE red
Diameter in mm	$8,30 \pm 0,15$
Total thickness in mm	$1,00 \pm 0,15$
Thickness PTFE	0,08
Durometer	45 ± 5

This document was sent online and is valid without signature