



Article Specification / Combination Seal

LLG-Screw caps N9, blue PP, closed top, Silicone white/PTFE red,

hardness: 40° shore A, thickness: 1.0 mm, pack of 100

Art.-Nr.: 4.008 227

Kappe/Cap

Kappentyp/Gewindetyp: Schraubkappe
Type of Cap/Type of Thread: Screw Cap

Material: Polypropylen
Material: Polypropylene

Legierung / Alloy

Lackierung/Farbe: Blau
Lacquer/Colour: blue

Abmessung/Dimensions: in mm

Kappenhöhe/Skirt Height: 5.85 ± 0,15

Materialdicke/Material Thickness:

Innendurchmesser/Inner Diameter: 8.70 ± 0,1

Außendurchmesser/Outer Diameter: 11.50 ± 0,1

Lochgröße/Size of Centre Hole:

Dichtscheibe/Septum

Material und Farbe: Silicon weiß/PTFE rot
Material and colour: Silicone white/PTFE red

Durchmesser/Diameter in mm: for pressfit

Gesamtdicke/Total thickness in mm: 1.00 ± 0,13

PTFE-Dicke/Thickness PTFE in mm: 0,08

Durometer/Durometer in shore A/D: 55 ± 5

This document was sent online and is valid without signature