



ARTICLE SPECIFICATION

LLG-Headspace Rollrandflasche N 20,
5ml A.D.: 21,7 mm, Außenhöhe: 38,2 mm, klar,
flacher Boden, abgeschrägter Rollrand, VE=100
Cat. No.: 4.008 285

Material:	Glass, 1st hydrolytic class
Colour:	Clear
Label:	No
Filling Lines:	No
Logo:	No
Dimensions [mm]	
Vial Height [mm]	38,25 ±0,5
Outer Diameter [mm]	21,75 ±0.25
Inner Diameter Bottleneck [mm]	12,60 ±0.2

This document was sent online and is valid without signature