



PRODUCT DATA SHEET

LLG-Centrifuge tubes economy

Cat. No.: 4.668 482, 4.668 483, 4.668 486, 6.270 403, 6.270 404 and 6.270 405



Description

15mL, Conical Bottom, Flat Cap, Sterile/Non-sterile

50mL, Conical Bottom/Self-standing, Flat Cap, Sterile/Non-sterile

Purpose

Popular in samples centrifuging or storage in many research area

Materials

Tube: PP (Polypropylene)

Color: Clear


Cap: HDPE (High-density Polyethylene)

Color: Orange

Features

- Autoclavable at 121°C and freezable to -80°C
- Easy-to-read black graduations in $\pm 2\%$ increments
- Contains a large, white frosted writing area
- Max rotate speed up to 12,000xg for conical bottom tubes, while 6,000xg for self-standing tubes
- DNase/RNase free
- Non-pyrogenic
- Sterilized by irradiation $SAL10^{-6}$ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system
- When the centrifuge tube is frozen, it is prohibited to use a foam shelf to freeze storage

Dimensions (Unit: mm)



Cat. No.	Capacity	Bottom Style	Cap Style	Sterilization	Packaging Configuration
4.668 482	15mL	Conical	Flat Cap	Yes	25/bag, 500/case
4.668 483	50mL	Conical	Flat Cap	Yes	25/bag, 500/case
4.668 486	50mL	Self-standing	Flat Cap	Yes	25/bag, 500/case
6.270 403	15mL	Conical	Flat Cap	No	500/bag/case
6.270 404	50mL	Conical	Flat Cap	No	500/bag/case
6.270 405	50mL	Self-standing	Flat Cap	No	500/bag/case

This document was sent online and is valid without signature