




Product Data Sheet

LLG-Crystallising dishes without spout, borosilicate glass 3.3

Cat. No.: 4.705 699



Article description: LLG-Crystallising Dishes without Spout, boro 3.3, 300 ml, 95x55 mm

The crystallising dishes have a flat base, are made of borosilicate glass 3.3, and are characterised by high temperature and chemical resistance. Complies with DIN 12337. Autoclavable.

Nominal capacity: 300 ml

Dimensions ($\varnothing$ x H): 95 x 55 mm

Packing unit: 1

This document was sent online and is valid without signature