




Product Data Sheet

LLG-Crystallising dishes without spout, borosilicate glass 3.3

Cat. No.: 4.705 701



Article description: LLG-Crystallising Dishes without Spout, boro 3.3, 900 ml, 140x75 mm

The crystallising dishes have a flat base, are made of borosilicate glass 3.3, and are characterised by high temperature and chemical resistance. Complies with DIN 12337. Autoclavable.

Nominal capacity: 900 ml

Dimensions (Ø x H): 140 x 75 mm

Packing unit: 1

This document was sent online and is valid without signature