



Acide chlorhydrique 1 mol/l (1N)

Synonyms:
HCl

Classification		Danger	
Molecular weight	ONU	1789	
CAS	7647-01-0	Transport hazard class:	8
EEC-N:	231-595-7	Packing	III

Acide chlorhydrique 1 mol/l (1N)

> RPE - Pour analyse

Grade Analytique

Description	Clear colourless liquid	-
Assay (potentiometry)	0.999 - 1.001	N
NIST 723	e	-
Description	Clear colourless liquid	-
Assay (potentiometry)	0.999 - 1.001	N
NIST 723	e	-



Acide chlorhydrique 1 mol/l (1N) > RPE - Pour analyse		Grade Analytique
ID	Size	Packaging
000000000404091000	10 l	Kubidos®
000000000404092000	5 l	Kubidos®
000000000404094000	20 l	Bidon plastique
000000000404097000	1l	Flacon plastique



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System