



Acide sulfurique 98%

Synonyms:

H₂SO₄**Classification****Danger**

Molecular weight

98.07 ONU**1830**

CAS

7664-93-9 Transport hazard
class:**8**

EEC-N:

231-639-5 Packing**II****Acide sulfurique 98%****> RS - Pour microanalyse****Grade Spécifique**

Description	Clear liquid	-
Identification	Positive	-
Density at 20°C	~ 1.835	-
Assay (acidimetric)	95 - 98	%
Description	Clear liquid	-
Identification	Positive	-
Density at 20°C	~ 1.835	-



Assay (acidimetric)	95 - 98	%
---------------------	---------	---

ID	Size	Packaging
000000000410421000	1 l	Flacon verre