



Butanol-2

Synonyms:

2-Butanol , alcool sec-butylique

CH₃CHOHCH₂CH₃

Classification

Danger

Molecular weight

74.12 ONU

1120



CAS

78-92-2 Transport hazard class:

3

EEC-N:

201-158-5 Packing

III

Butanol-2

> RPE - Pour analyse - Reag. Ph. Eur.

Grade Analytique

Description	Clear liquid	-
Colour	<= 10	APHA
Identification	Positive	-
Density d _{20/4}	0.801 - 0.811	-
Refractive index at 20°C	1.3944 - 1.3984	-
Boiling point	99 - 100	?C
Water (K.F.)	<= 0.2	%



Butanol-2

> RPE - Pour analyse - Reag. Ph. Eur.

Grade Analytique

Residue on evaporation	<= 20	ppm
Acidity (butyric acid)	<= 20	ppm
Alcalinity (NaOH)	<= 10	ppm
Indole base (indole)	<= 0.1	ppm
Assay (GC)	>= 99.0	%
Description	Clear liquid	-
Colour	<= 10	APHA
Identification	Positive	-
Density d20/4	0.801 - 0.811	-
Refractive index at 20°C	1.3944 - 1.3984	-
Boiling point	99 - 100	°C
Water (K.F.)	<= 0.2	%
Residue on evaporation	<= 20	ppm
Acidity (butyric acid)	<= 20	ppm
Alcalinity (NaOH)	<= 10	ppm
Indole base (indole)	<= 0.1	ppm
Assay (GC)	>= 99.0	%



ID	Size	Packaging
0000000000000414261	160 kg	Fût métal
0000000000000414264	1 l	Flacon verre
0000000000000414266	15 kg	Fût métal