



Ammoniaque solution 25%

Synonyms:

Classification		Danger	
Molecular weight	ONU	2672	
CAS	1336-21-6	Transport hazard class:	8
EEC-N:	215-647-6	Packing	III

Ammoniaque solution
25%

> RS - RSE - Pour usage électronique

Grade Spécifique

Description	Clear liquid	-
Colour	≤ 10	APHA
Identification	Positive	-
Density d20/4	0.904 - 0.910	-
Carbonate	≤ 5	ppm
Chloride	≤ 0.25	ppm
Phosphate	≤ 0.2	ppm



**Ammoniaque solution
25%****> RS - RSE - Pour usage électronique****Grade Spécifique**

Heavy metals (as Pb)	≤ 0.2	ppm
Residue on ignition	≤ 3	ppm
Subs. reducing KMnO ₄	≤ 5	ppm
Total sulphur (S)	≤ 1	ppm
Ag	≤ 0.02	ppm
Al	≤ 0.05	ppm
As	≤ 0.025	ppm
Au	≤ 0.05	ppm
B	≤ 0.01	ppm
Ba	≤ 0.1	ppm
Be	≤ 0.02	ppm
Bi	≤ 0.02	ppm
Ca	≤ 0.2	ppm
Cd	≤ 0.01	ppm
Co	≤ 0.01	ppm
Cr	≤ 0.01	ppm
Cu	≤ 0.01	ppm
Fe	≤ 0.03	ppm
Ga	≤ 0.02	ppm
In	≤ 0.02	ppm



K	<= 0.2	ppm
Li	<= 0.02	ppm
Mg	<= 0.1	ppm
Mn	<= 0.01	ppm
Mo	<= 0.05	ppm
Na	<= 0.5	ppm
Ni	<= 0.01	ppm
Pb	<= 0.01	ppm
Pt	<= 0.05	ppm
Sb	<= 0.05	ppm
Sn	<= 0.02	ppm
Sr	<= 0.02	ppm
Ta	<= 0.1	ppm
Ti	<= 0.05	ppm
Tl	<= 0.05	ppm
V	<= 0.05	ppm
Zn	<= 0.05	ppm
Zr	<= 0.05	ppm
Assay (alcalimetric)	24.0 - 26.0	%
Description	Clear liquid	-



**Ammoniaque solution
25%****> RS - RSE - Pour usage électronique****Grade Spécifique**

Colour	<= 10	APHA
Identification	Positive	-
Density d20/4	0.904 - 0.910	-
Carbonate	<= 5	ppm
Chloride	<= 0.25	ppm
Phosphate	<= 0.2	ppm
Heavy metals (as Pb)	<= 0.2	ppm
Residue on ignition	<= 3	ppm
Subs. reducing KMnO ₄	<= 5	ppm
Total sulphur (S)	<= 1	ppm
Ag	<= 0.02	ppm
Al	<= 0.05	ppm
As	<= 0.025	ppm
Au	<= 0.05	ppm
B	<= 0.01	ppm
Ba	<= 0.1	ppm
Be	<= 0.02	ppm
Bi	<= 0.02	ppm
Ca	<= 0.2	ppm
Cd	<= 0.01	ppm



**Ammoniaque solution
25%****> RS - RSE - Pour usage électronique****Grade Specific**

Co	≤ 0.01	ppm
Cr	≤ 0.01	ppm
Cu	≤ 0.01	ppm
Fe	≤ 0.03	ppm
Ga	≤ 0.02	ppm
In	≤ 0.02	ppm
K	≤ 0.2	ppm
Li	≤ 0.02	ppm
Mg	≤ 0.1	ppm
Mn	≤ 0.01	ppm
Mo	≤ 0.05	ppm
Na	≤ 0.5	ppm
Ni	≤ 0.01	ppm
Pb	≤ 0.01	ppm
Pt	≤ 0.05	ppm
Sb	≤ 0.05	ppm
Sn	≤ 0.02	ppm
Sr	≤ 0.02	ppm
Ta	≤ 0.1	ppm
Ti	≤ 0.05	ppm



Tl	<= 0.05	ppm
V	<= 0.05	ppm
Zn	<= 0.05	ppm
Zr	<= 0.05	ppm
Assay (alcalimetric)	24.0 - 26.0	%

ID	Size	Packaging
0000000000000420081	1 l	Flacon plastique
0000000000000420085	5 l	Flacon plastique



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System

