



Ether de pétrole 60 -80°C

Synonyms:

	Classification	Danger
Molecular weight	ONU	1268  
CAS	Transport hazard class:	3  
EEC-N:	Packing	II

Ether de pétrole 60 -80°C > RPE - Pour analyse

Grade Analytique

Description	Clear liquid	-
Colour	<=10	APHA
Identification	Positive	-
Boiling point	60 - 80	?C
Readily carbonizable substances	Conform	-
Density at 20°C	0.660 - 0.690	-
Water (K.F.)	<=100	ppm



Residue on evaporation	<=10	ppm
Acidity (acetic acid)	<=5	ppm
Subs. reducing KMnO4	<=20	ppm (5m)
Total sulphur	<=50	ppm
Description	Clear liquid	-
Colour	<=10	APHA
Identification	Positive	-
Boiling point	60 - 80	?C
Readily carbonizable substances	Conform	-
Density at 20°C	0.660 - 0.690	-
Water (K.F.)	<=100	ppm
Residue on evaporation	<=10	ppm
Acidity (acetic acid)	<=5	ppm
Subs. reducing KMnO4	<=20	ppm (5m)
Total sulphur	<=50	ppm



Ether de pétrole 60 -80°C > RPE - Pour analyse

Grade Analytique

ID	Size	Packaging
0000000000000427001	1 l	Flacon verre
0000000000000427003	2.5 l	Flacon verre
0000000000000427007	18 kg	Fût métal



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System

