



Ether de pétrole 80 - 100°C

Synonyms:

Classification		Danger	
Molecular weight	ONU	1268	 
CAS	64742-49-0	Transport hazard class:	3  
EEC-N:	Packing	II	

Ether de pétrole 80 - 100°C

> RPE - Pour analyse

Grade Analytique

Description	Clear colourless liquid	-
Colour	<=10	APHA
Identification	Positive	-
Boiling point	80 - 100	?C
Readily carbonizable substances	Conform	-
Density at 15°C	0.690 - 0.710	-
Water (K.F.)	<=0.01	%



**Ether de pétrole 80 -
100°C****> RPE - Pour analyse****Grade Analytique**

Residue on evaporation	<=10	ppm
Acidity (acetic acid)	<=5	ppm
Aromatic compounds	<=100	ppm
Subs. reducing KMnO ₄	<=20	ppm (5m)
Total sulphur	<=50	ppm
Description	Clear colourless liquid	-
Colour	<=10	APHA
Identification	Positive	-
Boiling point	80 - 100	?C
Readily carbonizable substances	Conform	-
Density at 15°C	0.690 - 0.710	-
Water (K.F.)	<=0.01	%
Residue on evaporation	<=10	ppm
Acidity (acetic acid)	<=5	ppm
Aromatic compounds	<=100	ppm
Subs. reducing KMnO ₄	<=20	ppm (5m)
Total sulphur	<=50	ppm



Ether de pétrole 80 - 100°C		> RPE - Pour analyse	Grade Analytique
ID		Size	Packaging
0000000000000427031		1 l	Flacon verre
0000000000000427036		18 kg	Fût métal