



Phénol

Synonyms:

Hydroxybenzène

C₆H₅OH**Classification****Danger**

Molecular weight

94.11 ONU**1671**

CAS

108-95-2 Transport hazard
class:**6.1**

EEC-N:

203-632-7 Packing**II****Phénol****> RPE - Pour analyse - ACS****Grade Analytique**

Description	Crystalline mass	-
Identification	Positive	-
Appearance of solution	Conform	ACS
Freezing point	>= 40.5	?C
Water (K.F.)	<= 0.5	%
Residue after evaporation	<= 0.05	%
Assay	>= 99.0	%



Description	Crystalline mass	-
Identification	Positive	-
Appearance of solution	Conform	ACS
Freezing point	>= 40.5	?C
Water (K.F.)	<= 0.5	%
Residue after evaporation	<= 0.05	%
Assay	>= 99.0	%

ID	Size	Packaging
000000000000451271	1 kg	Flacon plastique



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System

