



Jaune d'alizarine R

Synonyms:

Mordant Orange 1 , acide 5- (4-nitrophénylazo) salicylique

C13H9N3O5

Classification

Danger

Molecular weight **287.23** ONU



CAS **2243-76-7** Transport hazard class:

EEC-N: **218-818-3** Packing

Jaune d'alizarine R

> RPE - Pour analyse - C.I. 14030

Grade Analytique

Description	Brown crystalline powder	-
Identification	Positive	-
Loss on drying	<= 15	%
Colour change	Yellow-orange	-
E (1% 1cm) at 492nm	700 - 1000	-



ID	Size	Packaging
000000000000453451	10 g	Flacon verre
000000000000453452	25 g	Flacon verre