



8-Hydroxyquinoléine

Synonyms:

HOC6H3N:CHCH:CH**Classification****Danger**

Molecular weight

145.16 ONU**2811**

CAS

148-24-3 Transport hazard
class:**6.1**

EEC-N:

205-711-1 Packing**III****8-Hydroxyquinoléine****> RPE - Pour analyse****Grade Analytique**

Description	Off white to light brown crystal	-
Identification	Positive	-
Melting point	72.5 - 74.5	?C
Residue on ignition	<= 2000	ppm
Sulfate	<= 200	ppm
Assay	>= 99.5	%



Description	Off white to light brown crystal	-
Identification	Positive	-
Melting point	72.5 - 74.5	?C
Residue on ignition	<= 2000	ppm
Sulfate	<= 200	ppm
Assay	>= 99.5	%

ID	Size	Packaging
0000000000000467353	50 g	Flacon verre
0000000000000467355	250 g	Flacon plastique
0000000000000467356	1 kg	Flacon plastique