



Acide chlorhydrique 6 mol/l (6N)

Synonyms:

HCl

Classification

Danger

Molecular weight

ONU

1789



CAS

7647-01-0 Transport hazard
class:

8

EEC-N:

231-595-7 Packing

III

**Acide chlorhydrique 6
mol/l (6N)**

> RS - Pour analyse agroalimentaire

Grade Spécifique

Description	Clear colourless liquid	-
Assay (potentiometry)	5.97 - 6.03	N
Colour	<= 10	APHA
Assay (potentiometry)	5.97 - 6.03	N
Colour	<= 10	APHA
Description	Clear colourless liquid	-



Acide chlorhydrique 6 mol/l (6N) > RS - Pour analyse agroalimentaire		Grade Specifique
ID	Size	Packaging
000000000000502831	1 l	Flacon plastique
000000000000502832	18 l	Bidon plastique