



Acide chlorhydrique 50% v/v

Synonyms:

HCl

Classification

Danger

Molecular weight	MM	ONU	1789		
	36.461				
CAS	7647-01-0	Transport hazard class:	8		
EEC-N:	231-595-7	Packing	III		

Acide chlorhydrique 50% v/v > RPE - Pour analyse

Grade Analytique

Density d20/4	1.092 - 1.102	-
Density d20/4	1.092 - 1.102	-



Acide chlorhydrique 50% > RPE - Pour analyse
v/v

Grade Analytique

ID	Size	Packaging
000000000000504571	1 l	Flacon verre



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System

