



Propane-2-ol

Synonyms:

2-Propanol , Isopropanol , IPA

CH3CHOHCH3

Classification

Danger

Molecular weight

60 ONU

1219



CAS

67-63-0 Transport hazard
class:

3

EEC-N:

200-661-7 Packing

II

Propane-2-ol

> RS - VLSI Pour usage électronique

Grade Specificque

Assay	>= 99.9	%
Density at 20°C	~ 0.785	-
Colour	<= 10	APHA
Dry residu	<= 5	ppm
Chloride	<= 1	ppm
Nitrate	<= 1	ppm
Sulfate	<= 1	ppm



Sulfite	<= 1	ppm
Acidity	<= 0.2	?eq/g
Alcalinity	<= 0.1	?eq/g
Water (K.F.)	<= 1000	ppm
Other impurities	-	-
Al	<= 20	ppb
As	<= 10	ppb
Bi	<= 10	ppb
Cd	<= 10	ppb
Cr	<= 10	ppb
Cu	<= 10	ppb
Au	<= 20	ppb
Li	<= 10	ppb
Mn	<= 10	ppb
Ni	<= 10	ppb
K	<= 20	ppb
Na	<= 50	ppb
Sn	<= 10	ppb
V	<= 10	ppb
Sb	<= 10	ppb



Ba	<= 20	ppb
B	<= 20	ppb
Ca	<= 50	ppb
Co	<= 10	ppb
Ga	<= 10	ppb
Fe	<= 50	ppb
Mg	<= 50	ppb
Mo	<= 10	ppb
Pb	<= 10	ppb
Ag	<= 10	ppb
Sr	<= 20	ppb
Ti	<= 10	ppb
Zn	<= 20	ppb
Assay	>= 99.9	%
Density at 20°C	~ 0.785	-
Colour	<= 10	APHA
Dry residu	<= 5	ppm
Chloride	<= 1	ppm
Nitrate	<= 1	ppm
Sulfate	<= 1	ppm



Sulfite	<= 1	ppm
Acidity	<= 0.2	?eq/g
Alcalinity	<= 0.1	?eq/g
Water (K.F.)	<= 1000	ppm
Other impurities	-	-
Al	<= 20	ppb
As	<= 10	ppb
Bi	<= 10	ppb
Cd	<= 10	ppb
Cr	<= 10	ppb
Cu	<= 10	ppb
Au	<= 20	ppb
Li	<= 10	ppb
Mn	<= 10	ppb
Ni	<= 10	ppb
K	<= 20	ppb
Na	<= 50	ppb
Sn	<= 10	ppb
V	<= 10	ppb
Sb	<= 10	ppb



Ba	<= 20	ppb
B	<= 20	ppb
Ca	<= 50	ppb
Co	<= 10	ppb
Ga	<= 10	ppb
Fe	<= 50	ppb
Mg	<= 50	ppb
Mo	<= 10	ppb
Pb	<= 10	ppb
Ag	<= 10	ppb
Sr	<= 20	ppb
Ti	<= 10	ppb
Zn	<= 20	ppb

ID	Size	Packaging
0000000000000527690	2.5 l	Flacon plastique
0000000000000527693	25 l	Bidon plastique
0000000000000527696	1 l	Flacon plastique