



Ethanol 70% v/v

Synonyms:

	Classification	Danger
Molecular weight	ONU	1170  
CAS	64-17-5 Transport hazard class:	3
EEC-N:	Packing	II

Ethanol 70% v/v

> RE - Puro

Grade Technique

Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Residue on evaporation	<= 50	ppm
Assay (alcoholic) at 20°C	69 - 71	% (V/V)
Residue on evaporation	<= 50	ppm
Assay (alcoholic) at 20°C	69 - 71	% (V/V)
Description	Clear colourless liquid	-



Identification (I.R.)	Positive	-
-----------------------	----------	---

ID	Size	Packaging
000000000000528170	5 l	Bidon plastique
000000000000528171	1 l	Flacon plastique
000000000000528172	2.5 l	flacon plastique
000000000000528173	10 l	bidon plastique
000000000000528174	25 l	bidon plastique