



Ammoniaque solution 32%

Synonyms:

Classification		Danger	
Molecular weight	ONU	2672	 
CAS	1336-21-6	Transport hazard class:	8 
EEC-N:	215-647-6	Packing	III

**Ammoniaque solution
32%**

> RE - Puro

Grade Technique

Description	Clear liquid	-
Colour	<= 10	APHA
Identification	Positive	-
Density at 20°C	0.880 - 0.898	-
Assay (alcalimetric)	28 - 34	%
Description	Clear liquid	-
Colour	<= 10	APHA



Ammoniaque solution > RE - Puro		Grade Technique
32%		
Identification	Positive	-
Density at 20°C	0.880 - 0.898	-
Assay (alcalimetric)	28 - 34	%

ID	Size	Packaging
000000000000528501	18 kg	Fût combiné métal / plastique
000000000000528503	5 l	Flacon plastique