



Ethanol 70% v/v

Synonyms:

Classification

Danger

Molecular weight **46.070** ONU**1170**  CAS **64-17-5** Transport hazard class:**LQ**

EEC-N: Packing

Ethanol 70% v/v

> ERBApharm - Selon pharmacopée: Ph.Eur.

Grade Pharmaceutique

Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Assay (alcoholic) at 20°C	69 - 71	% (V/V)
Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Assay (alcoholic) at 20°C	69 - 71	% (V/V)



Ethanol 70% v/v

> ERBApharm - Selon pharmacopée: Ph.Eur.

Grade Pharmaceutique

ID	Size	Packaging
000000000000529183	200 l	Fût plastique
000000000000529184	1 l	Flacon spray
000000000000529189	5 l	Bidon plastique
000000000529184000	1 l	Flacon spray
000000000529187000	500 ml	Flacon spray



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val de
Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA Reagents
operates with a Certified
Quality Management
System