

## CERTIFICATE OF QUALITY CONTROL

**Description** : Measuring Cylinders - "embossed graduation"  
**Catalog Number** : 016.05.901  
**Material** : Polypropylene  
**Color** : Clear  
**Dimensional Control** : Pass  
**Volumetric Control** : Pass  
**Visual Appearance** : Within acceptability limits  
**Calibration** : Class B tolerances according to DIN 12681 and ISO 6706

Catalog Number	Volume	Subdivision	Tolerance	Approx. Height (h)
	ml	ml	ml	mm
016.05.901	1000	10.000	± 10.000	465

We declare that the above mentioned batch on this certificate is manufactured under the ISO 9001:2008 Quality Management System. Products meet the highest quality standards according to DIN 12681 and ISO 6706. Cylinders are calibrated to contain (In) according to Class B tolerances. Measuring cylinders should not be kept under direct sunlight or overheated conditions.

The above mentioned properties tested by ISOLAB Quality Control Laboratory.

For and on behalf of

**Quality Assurance Department**  
**ISOLAB LABORGERÄTE GmbH**

**Date** : August 2015

**Quality Assurance Manager**

