



ARTICLE SPECIFICATION

LLG-Magnetic screw cap N 18, silver, center hole Silicone blue transparent/PTFE white,
Hardness: 45° shore A, Thickness: 1.5 mm

Cat. No.: 6.241 111

Cap	
Type of Cap/Typ of Thread:	magnetic Universal-Screw Cap
Lacquer/Colour:	silver
Dimensions mm:	
Skirt Height :	11.00 ± 0.3
Material Thickness :	0.20 ± 0.01
Inner Diameter :	16.80 + 0.17/-0.1
Outer Diameter :	18.40 ± 0.12
Size of Centre Hole :	8.00 ± 0.1
Septum	
Material and colour:	Silicone transparent blue / PTFE white
Diameter in mm:	for pressfit
Total thickness in mm:	1.30 ± 0.15
Thickness PTFE in mm:	0.13
Durometer in shore A/D:	45 ± 5

This document was sent online and is valid without signature