



Technical Data Sheet

LLG-Thread bottles 8 ml thread 15-425, pack of 100

Art.-Nr.: 6.280 953

Type of Cap: 8 ml Screw Neck Vial, 15-425 thread

Material: Glass, 1st hydrolytic class

Colour: Clear

Label: No

Filling Lines: No

Logo: No

Dimensions [mm]

Vial Height: 61.00 $\pm$ 0.5

Outer Diameter: 16.60 $\pm$ 0.25

Inner Diameter Bottleneck: 13.05 $\pm$ 0.20

This document was sent online and is valid without signature