



TECHNICAL DATA SHEET

LLG-short thread caps 9 mm PP, blue, with hole, silicone white/ PTFE red, pack of 100

Cat. No.: 6.291 638

Cap

Type of Cap: Screw Cap with Centre Hole
Material: Polypropylene
Lacquer: blue

Dimensions:

Skirt Height: $5,70 \pm 0,20$
Material Thickness: $0,27 \pm 0,02$
Inner Diameter: $11,00 \pm 0,10$
Outer Diameter: $11,5 \pm 0,12$
Size of Centre Hole: $5,55$

Septum

Material: Silicone white/PTFE red, 45° shore A, 1.0 mm
Colour : white/red
Diameter : $9,0 \pm 0,15$ mm
Total Thickn. : $1,0 \pm 0,15$ mm
Thickn. PTFE : 0,1 mm
Durometer : 45 $\pm$ 5° shore A
Special Properties: none

This document was sent online and is valid without signature