



Product data sheet

LLG-Protection glasses "Comfort"

Art.-N.: 6. 312 450



Frame mark (EU)

U EN166 FT CE

Lens mark (EU)

2C-1.2 U 1 FT K N CE

Frame technology

Adjustment of temples in length, Adjustment of temples in inclination, Rubber Terminal, Terminal tilting movement, SoftPad

Conditions of use

Standard brightness,
Low brightness

Lens technology

Anti-Scratch + Anti-Fog
(KN - Vanguard PLUS), UV400

Lens color

Clear

Lens material

Polycarbonate

Frame color

Gunmetal/Green

Standards

EN166, EN170

Scale number

(2C-1.2)

Mechanical resistance

Low energy impact (F)

Optional requirements

Resistance to surface damage by fine particles (K), Resistance to fogging of oculars (N), Protection against high speed particles at extreme temperatures -5°/+55° C (T)

Optical class

Optical class (1)

Weight

30 g

Reference standards: EN 166, EN 170

This document was sent online and is valid without signature