



DURAN®

AUTOCLAVE BOTTLE

The DURAN® Autoclave Bottle is the preferred method for sterilizing water without the need for post-autoclaving intervention. Additionally, when paired with the DURAN® 45mm self-venting natural rubber cap, ongoing sterility integrity is visually assured thanks to the physical shape of the concave cap.

PRODUCT ADVANTAGES

- Fully compatible with 45mm self-venting, natural rubber sealing cap DURAN® closure, the preferred solution for hands-free sterilization of lab water
- Autoclavable at 121°C or 134 °C
- High-purity, clear, USP/EP/JP Type I glass
- Borosilicate glass 3.3, renowned for its excellent thermal performance, and high resistance to chemical attack
- Retrace code for manufacturing lot traceability
- Fired-on, highly durable white print
- Easy-to-read double graduated scales to show either liquid volume contained or volume poured out
- Single-bottle packaged for easy ordering, handling, and inventory

ORDERING INFORMATION

| Cat. No. | Description | Qty / Case |
|-----------|--|------------|
| 219946305 | DURAN® Autoclave Bottle 2000mL | 1 |
| 292063002 | DURAN® 45mm Push-on Natural Rubber Cap | 100 |

