



ARTICLE SPECIFICATION

LLG-Probengläser N 20-6 PE, farblos Rollrand, VE=100

Cat. No.: 7.052 186

Material:	Glass, 1st hydrolytic class
Colour:	Clear
Label:	No
Filling Lines:	No
Logo:	No
Dimensions [mm]	
Vial Height [mm]	38,25 ±0,5
Outer Diameter [mm]	22,00 ±0.25
Inner Diameter Bottleneck [mm]	12,60 ±0.20

This document was sent online and is valid without signature