



ARTICLE SPECIFICATION

LLG-Screw Neck Vial N 13, 4ml, O.D.: 14.75mm outer height: 45 mm, amber,
flat bottom pack of 100
Cat. No.: 7.058 142

Material:	Glas, 1 st hydrolytic class
Colour:	Amber
Label:	No
Filling Lines:	No
Logo:	No
Dimensions [mm]	
Vial Height:	45.00 ± 0.5
Outer Diameter:	14.75 ± 0.15
Inner Diameter Bottleneck:	8.50 ± 0.15

This document was sent online and is valid without signature