



TECHNICAL DATA SHEET

LLG-Headspace-Vial ND20, 10 ml, 20.0 x 54.5 mm, Brown, DIN crimp neck, flat, Varian

pack of 100pcs

Art.-Nr.: 7.080 947



Parameter*	Value*
Silanized	No
Usable volume	10.0 mL
Vial category	Vial
Vial diameter	20.5 mm
Vial dimensions	20.5 mm x 54.5 mm
Vial height	54.5 mm
Vial material	Amber glass, 1st hydrolytic class glass
Bottle type	Crimp neck vial
Bottom shape of the vial	Flat bottom
Writing patch	No
Additional vial features	Flat DIN crimp neck
Autosampler compatibility	Varian 8410 Autoinjector
Nominal diameter	N 20
Recommended application(s)	HPLC, Headspace
Scope of delivery	100x N 20 glass crimp neck vials
Storage temperature	15–25 °C / 59–77 °F
Hazardous material	No Information available
Gross weight (incl. packaging)	1027.5 g
Packaging dimensions	190 x 245 x 60 mm

*These data are typical values and must not be regarded as specification

This document was sent online and is valid without signature