



ARTICLE SPECIFICATION

LLG-Screw cap N 13, PP, black, center hole, Silicone, white/PTFE
red, Hardness: 40° shore A, Thickness: 1.5 mm
Cat. No.: 7.510 053

Cap	
Material:	Polypropylene
Colour:	black
Dimensions [mm]	
Skirt Height:	9,45 ± 0,15
Inner Diameter:	11,90 ± 0,15
Outer Diameter:	16,25 ± 0,2
Size of Centre Hole :	8,70 ± 0,2
Septum	
Material and colour:	Silicone cream / PTFE red
Diameter in mm:	for pressfit
Total thickness in mm:	1,50 ± 0,13
Thickness PTFE in mm:	0,08 +0,01
Durometer in shore A/D:	55 ± 5

This document was sent online and is valid without signature