



SPECIFICATION SHEET

Art.-Nr.: 7.606 811

DESCRIPTION


Divider made of paperboard with 5 x 5 cells.
Designed for cryo vials smaller than 24 mm.
Divider is 30 mm high and fits to a cryobox 50 mm high.
Designed for long-term storage of samples at ultra low temperatures.
Made in EU.

MATERIALS

PAPERBOARD

- double coated paperboard
- basis weight is 300 g / m² +-4% (ISO536)
- thickness is 0,474 mm +-4% (ISO534)
- moisture content 7,8% +-1,0% (ISO287)
- all fiber originates from sustainably managed forests
- produced at a mill that is certified to the ISO14001 environmental management standard
- FSC mix credit system (SGSCOCO02560)
- ISEGA certificate allowing food applications
- PTS / FFPI certification for inkjet coding on pharmaceutical packaging.

TECHNICAL SPECIFICATIONS



EXTERNAL DIMENSIONS:
123 x 123 x 30 mm

CELL DIMENSIONS:
24,6 x 24,6 mm

NET WEIGHT OF THE PRODUCT:
8 g

PACKAGE: 36 pcs / carton

CARTON EXTERNAL DIMENSIONS:
250 x 240 x 78 mm

GROSS WEIGHT: 409 g
NET WEIGHT: 288 g

Lab Logistics Group GmbH
Product Management
Am Hambuch 1
D-53340 Meckenheim / Germany

Meckenheim, 2014-03-27

This document was sent online and is valid without signature

