



Technical Data Sheet

Product No: 7.930 366

Product description: LLG Labware - 9mm Combination Seal: PP Short Thread Cap, blue, closed top. PTFE virginal, 53° shore D, 0,2 mm, pack of 100pcs

Type of Cap:	Screw Cap with Closed Top		
Material:	Polypropylene		
Lacquer:	blue		
Dimensions:	Skirt Height:	5,70 ± 0,20	
	Material Thickness:	0,27 ± 0,02	
	Inner Diameter:	11,0 ± 0,10	
	Outer Diameter:	11,5 ± 0,12	
	Size of Centre Hole:	5,55	

Product No: 4.714 105

Material:	9mm Septa, PTFE virginal, 53° shore D, white, 0.20mm
Diameter:	9,0 ± 0,15 mm
Total Thickn.:	0,20 mm
Thickn. PTFE:	0,20 mm
Durometer:	53 ± 5° shore D
Special Properties:	none

This document was sent online and is valid without signature