



ARTICLE SPECIFICATION

LLG-Probengläser N 20-5 DIN, farblos farblos, Rollrand-Höhe 3 mm,
5 ml 20,5 x 38 mm, VE=100
Cat. No.: 9.003 425

Material:	Glass, 1st hydrolytic class
Colour:	Clear
Label:	No
Filling Lines:	No
Logo:	No
Dimensions [mm]	
Vial Height [mm]	38,00 ±0,5
Outer Diameter [mm]	20,00 ±0.25
Inner Diameter Bottleneck [mm]	12,60 ±0.20

This document was sent online and is valid without signature