



TECHNICAL DATA SHEET

LLG-Crimp Seals ND20, Aluminium, ready assembled, silver, 10 mm centre hole, Material septa: Silicone/PTFE, Colour septa: Blue transparent/transparent, Shore hardness 45 A

pack of 100pcs

Art.-Nr.: 9.003 434

Type of Cap:	Crimp Cap
Material:	Aluminium
Colour:	clear lacquered
Skirt Height:	7.50 ± 0.15 mm
Material Thickness:	0.20 ± 0.01 mm
Inner Diameter:	20.30 ± 0.05 mm
Outer Diameter:	20.70 ± 0.07 mm
Size of Centre Hole:	10.00 ± 0.2 mm

This document was sent online and is valid without signature