



TECHNICAL PRODUCT SPECIFICATION

Cat. No.: 9006370 :: LLG-Disposable gloves, ERGO, Nitrile, size XS

Cat. No.: 9006371 :: LLG-Disposable gloves, ERGO, Nitrile, size S

Cat. No.: 9006372 :: LLG-Disposable gloves, ERGO, Nitrile, size M

Cat. No.: 9006373 :: LLG-Disposable gloves, ERGO, Nitrile, size L

Cat. No.: 9006374 :: LLG-Disposable gloves, ERGO, Nitrile, size XL

General information

Type	single use examination and disposable protective glove, non sterile
Labelling	information printed on dispenser box
Shape	ambidextrous - straight fingers
Material	Nitrile Butadiene Rubber (NBR) [not made with natural rubber latex]
Colour	violet blue
Inside	powder free
Outside	no treatment
Cuff / surface	rolled cuff / finger textured
Shelf life	5 years
Available sizes	XS (5-6) S (6-7) M (7-8) L (8-9) XL (9-10) XXL (10-11)

Dimensions, physical properties and biocompatibility

Glove length	median ≥ 240 mm	
Minimum wall thickness	<i>at finger</i>	0.14 mm (double measured) / 0.07 mm (single measured)
	<i>at palm</i>	0.12 mm (double measured) / 0.06 mm (single measured)
	<i>at cuff</i>	0.10 mm (double measured) / 0.05 mm (single measured)
Glove width	median XS ≤ 80 mm, S 80 ± 10 mm, M 95 ± 10 mm, L 110 ± 10 mm, XL ≥ 110 mm, XXL >115 mm	
Force at Break	median ≥ 6 N (during shelf life)	
Tensile Strength	min. 14 MPa after aging (according to ASTM D6319)	
Elongation at Break	min. 400% after aging (according to ASTM D6319)	
Residual powder / Powder content	≤ 2 mg	



Performance requirement and inspection level

Freedom from holes (Barrier)	AQL $\leq$ 1.0 sampling in accordance with ISO 2859-1, G-1)
-------------------------------------	--

Standards, guidelines & quality certificates

Quality certification	ISO 9001, ISO 14001
Conformity to regulations	- PPE Regulation (EU) 2016/425: Category III - Regulation (EC) 1935/2004 on Food Contact Materials
Conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, ISO 374-5, ISO 2859-1, ASTM F1671/F1671M

Instructions and additional statements

Storage instruction	Store in original packaging in a dry and dark place at 10° C to 30° C. Refer to guidelines of storage of rubber products as described in ISO 2230:2002. Ensure that storage area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper ions discolour the glove. Protect gloves against ultraviolet light sources, such as sunlight and oxidizing agents. Storage above 30° C will lead to accelerated aging and should be avoided.
Cautionary statement and ingredient information	Contains potential Type IV chemical allergens (Dithiocarbamate types). For further information, a list of substances contained in the glove is available upon request.

Date of issue: 08.10.2024

This document was sent online and is valid without signature

Page 2 of 2