




TECHNICAL DATA SHEET

LLG-Porcelain crucibles, tal form

The glazed porcelain crucible is characterised by high thermal and chemical resistance. Ideal for heating and melting of materials. According to DIN 12904. Temperature resistant up to 900 °C



Capacity ml	Ext. Diam. mm	Height mm	Cat. No.
35	40	50	6.238 374
15	30	38	9.250 920
26	35	44	9.250 921
35	40	50	9.250 922
50	45	56	9.250 923
72	50	62	9.250 924
130	60	75	9.250 925
15	30	38	9.250 970
26	35	44	9.250 971

This document was sent online and is valid without signature