



Product data sheet

LLG-Flat bags, PE, 150 x 200 mm, pack of 100 Art.-Nr.: 9.404 130

Technical data:

Designation	Polybag unprinted
Opening	narrow side
Dimensions	according to delivery documents
Material thickness	50 my
Dimensions/Material thickness	Production is carried out according to the guidelines of the GKV
Material structure	LDPE, extruded, unstretched quality
Color	transparent
Closure	without
Slide description	Low-density Polyethylene
Country of origin	batch dependent on request
Customs tariff number	39232100000
Motif/Print	without
Test method	DIN53370
execution	unprinted
Sealability	thermo-weldable
Physiological safety	food safe according to BgVV, EC guidelines
Storage	cool and dry at 15-25 °C, 50% relative humidity (the original packaging offers the best protection)

Notes:

These technical specifications are intended as guidelines based on our knowledge and experience. They do not represent a guarantee of specific properties.

The packaging company is responsible for checking whether the packaging is suitable for the intended purpose. This applies in particular to changes in the quality of the contents that may arise from interaction with the packaging material or parts of it.

This document was sent online and is valid without signature